

PsychoPy:反应盒精确计时/Cedrus Box

Jibo He, Ph.D.

Department of Psychology, Tsinghua University

hejibolaboratory@pku.org.cn

反应盒与键盘Response box vs. keyboard

两者的反应时精度的差异

2-3ms vs. 10-12ms

" Get 2-3 millisecond RT resolution. The resolution of a USB keyboard driver is 10-12ms and the hardware could make it worse."

https://www.cedrus.com/rb_series/index.htm

反应盒与键盘Response box vs. keyboard

Table 2 Descriptive statistics for RTs and accuracy, aggregated by participants												
Set Size	Response Box			Keyboard			Web From Lab			Web From Home		
	2	4	6	2	4	6	2	4	6	2	4	6
RT mean	319	473	535	378	524	565	443	577	614	429	545	608
RT <i>SD</i>	66	114	78	60	94	46	133	215	105	72	94	87
Share of errors	4.3	6.0	6.4	2.5	3.8	5.1	3.7	4.5	5.4	4.0	5.2	6.2

Source: Plant, R. R. (2016).

反应盒与键盘Response box vs. keyboard

Behav Res (2016) 48:1086–1099

1091




Fig.2 Means (left panel) and distributions (right panel) of response times (RTs, in ms). Lines in the left panel show 95 % confidence intervals, controlling for between-subjects variability (Cousineau, 2005; Morey, 2008)


Source: Plant, R. R. (2016).

反应盒的类型


- Cedrus box (supported by PsychoPy)
- SRBox (for Eprime)

Cedrus 反应盒

- Cedrus 7-button response pad



Cedrus 反应盒



Cedrus 安装驱动

INSTALLING AN RB SERIES PAD

The RB-x40 generation response pads introduce a new Keyboard mode. For better timing, experimental software like SuperLab use Standard mode instead and need to have a USB driver installed:

- Installing the USB Driver for Mac OS X
- Installing the USB Driver for Windows
- Installing the USB Driver for Linux

<https://ftdichip.com/drivers/d2xx-drivers/>

Cedrus 提供的PsychoPy示例代码

https://cedrus.com/support/rb_series/tn1551_psychopy.htm

```
import pyxid2
import time

# get a list of all attached XID devices
devices = pyxid2.get_xid_devices()

dev = devices[0] # get the first device to use
print(dev)
dev.reset_base_timer()
dev.reset_rt_timer()

#dev._send_command('iuA1', 0)
#dev.enable_usb_output('K', True)

if dev.is_response_device():
    print("Press a key!")
    while not dev.has_response():
        dev.poll_for_response()


    response = dev.get_next_response()
    print(response)
    dev.clear_response_queue()

dev.set_pulse_duration(300)

sleep_flash = .3
for bm in range(0, 16):
    mask = 2 ** bm
    print("activate_line bitmask: ", mask)
    #dev.activate_line(lines=[1,3,5,7,9,11,13,15])
    dev.activate_line(bitmask=mask)

time.sleep(sleep_flash)
```

Run Cedrus box for PSYCHOPY in WINDOWS



The screenshot shows the PsychoPy Code Editor interface. The title bar indicates the file path: `\\Mac\Home\Downloads\psychopy_pyxid_test.py - PsychoPy Code`. The menu bar includes `File`, `Edit`, `View`, `Tools`, `Demos`, `Pavlovia.org`, `Window`, and `Help`. The toolbar contains icons for file operations and a green play button labeled `Run experiment`. The `Source Assistant` pane on the left shows the `psychopy_pyxid_test.py` file. The `Editor` pane displays the following Python code:

```
1 import pyxid2
2 import time
3
4 #get a list of all attached XID devices
5 devices = pyxid2.get_xid_devices()
6
7 dev = devices[0] #get the first device to use
8 print(dev)
9 dev.reset_base_timer()
10 dev.reset_rt_timer()
11
12 #dev._send_command('iuA1', 0)
13 #dev.enable_usb_output('K', True)
14
15 if dev.is_response_device():
16     print("Press a key!")
17 while not dev.has_response():
18     dev.poll_for_response()
```

The `Shelf` pane at the bottom shows the `PyShell` in `PsychoPy` with the prompt `>>>`. The status bar at the bottom indicates `Line: 1 Col: 1` and `Python`.

示例代码得到Cedrus box反应时

```
import pyxid2 as pyxid

# get a list of all attached XID devices
devices = pyxid.get_xid_devices()

dev = devices[0] # get the first device to use
if dev.is_response_device():
    dev.reset_base_timer()
    dev.reset_rt_timer()

    while True:
        dev.poll_for_response()
        if dev.response_queue_size() > 0:
            response = dev.get_next_response()
            # do something with the response
```

How to Cite Cedrus Box in your paper?

"The experiment was run on a Windows 7 PC. In the odor identification and rating tasks, participants responded with the mouse. In the categorization task, they responded with a Cedrus RB-740 Response Box (Cedrus Corporation)."

Source: Thomas Hörberg, et al.(2020)

使用反应盒的参考文献:

Thomas Hörberg, Maria Larsson, Ingrid Ekström, Camilla Sandöy, Peter Lundén, Jonas K Olofsson, Olfactory Influences on Visual Categorization: Behavioral and ERP Evidence, Cerebral Cortex, Volume 30, Issue 7, July 2020, Pages 4220–4237,
<https://doi.org/10.1093/cercor/bhaa050>

精确计时参考文献

Plant, R. R. (2016). A reminder on millisecond timing accuracy and potential replication failure in computer-based psychology experiments: An open letter. *Behavior Research Methods*, 48(1), 408-411.

References:

PsychoPy

- <https://www.psychopy.org/api/hardware/cedrus.html>

Cedrus

- <https://github.com/cedrus-opensource/pyxid>
- https://cedrus.com/support/rb_series/index.htm
- https://cedrus.com/support/rb_series/tn1551_psychopy.htm
- <https://phon.wordpress.ncsu.edu/lab-manual/perception-experiments/>